

AppBlaster User Guide

Version 1.0

11 May 2012

1

© CECID

Produced by the Center for E-Commerce Infrastructure Development
The University of Hong Kong

The contents of this document remain the property of and may not be reproduced in whole or in part without the express permission of CECID

1. Status of This Document

1.1. Version History


Ver. No.	Date	Revised By	Description	Filename

2. Table of Contents

1. STATUS OF THIS DOCUMENT	2
1.1. VERSION HISTORY	2
2. TABLE OF CONTENTS	3
3. INTRODUCTION	4
4. INSTALLATION & SETUP.....	5
4.1. SOFTWARE REQUIREMENT.....	5
4.2. INSTALLATION	5
4.3. SETUP	5
4.3.1. Bluetooth Connection setup with Control Gateway	5
4.3.2. Setup electric devices	5
4.3.3. Setup the camera zoom ratio	7
5. TRAINING PHASE.....	9
5.1. SETTING THE GAZE TIME FOR TRAINING.....	9
5.2. TRAINING.....	10
6. OPERATION PHASE.....	12
6.1. MOVE AND SELECTION.....	12
6.2. SET THE GAZE DURATION	13
6.3. LOCK SCREEN	13

3. Introduction

The AppBlaster is an android application that can enable disabilities to control the electric appliances through eye gaze. The following diagram illustrates how it works:



The App Blaster will send signal to Control Gateway using Bluetooth and then the Control Gateway will send the signal to the corresponding infrared controllable plugs, so that it can control the on/off for the electric appliances that plugged to them..

As the AppBlaster has to control the electric appliance, setup is necessary for configure the Bluetooth connection, and also what signal should be used to control electric appliances.

Please see the Installation & Setup section for details.

After the setup, users should first perform training phase on AppBlaster before the operation.

Please see the Training Phase section for details.

For the operation on AppBlaster, please see the Operation Phase for details.

4. Installation & Setup

4.1. Software Requirement

AppBlaster must be installed on the Android device with Android 2.2 or above. In order to have a better gaze direction recognition result, android tablet (around 10 inches) is preferred.

4.2. Installation

AppBlaster is an android application. Therefore the installation is easy, just like other android application.

4.3. Setup

4.3.1. Bluetooth Connection setup with Control Gateway

The AppBlaster is connected to Control Gateway through Bluetooth. In order to make the connection, it is necessary for the Android device to paired up and connect to the Control Gateway, which can be done in the setup on the Android device.

4.3.2. Setup electric devices

In this section, it will describe how to setup what signal should be used to control electric appliances.


When we start the AppBlaster, it will show the following screen:




Then user can click the setup menu button, the menu item will show.




- 1
- 2 Then user can click the setting button. The following screen will show.




- 3
- 4 Then user can click the “electric devices”, it will show the following screen.




- 5
- 6 The user can click the “add electric devices” button, it will show the screen to ask user for the
- 7 name of the electric device.



- 1
- 2 After the name is entered, user can enter the type and signal on the following screen.



- 3
- 4 The type means the type of the electric appliance. The signal means the signal that send to
- 5 Control Gateway. After setting the type and the signal, user can click ok to confirm the adding
- 6 of device.




- 7
- 8 **4.3.3. Setup the camera zoom ratio**

9 As for different device, the camera zoom ratio may be different. Therefore, user may choose

10 to change the camera zoom ratio for better eye detection and gaze direction detection


11 accuracy.

12 In order to change the camera zoom ratio, first go to the settings screen



1

2 Then user can select the camera zoom ratio to change the camera ratio.



3


4

5. Training Phase


In order to use the AppBlaster, user must at least perform the training phase once. In the training phase, AppBlaster will collect the eye gazing information from the user by asking the user to gaze on both the left and the right upper corner on the screen. The time required for gaze on each direction during the training can be configured. The longer the time, the more accurate for the training, however it may be rough to user.

5.1. Setting the gaze time for training


In order to set the time required on gazing each direction during training, the user can open the settings screen



Then the user click the training, the following screen will be shown.



Then the use click the gaze duration, it will allow user to choose the level for gaze duration.



- 1
2 After user selection, user may then exit the setting screen.

3 5.2. Training

- 4 When the AppBlaster starts, it will on the following screen




- 5
6 In the screen, there is a camera view. In the camera view there is a white rectangle. User
7 should try to move the head so that his/her eyes are within the center of the white rectangle. If
8 the system can detect the eyes, green rectangle will be shown. User should try to configure the
9 camera zoom ratio so that the green rectangle can enclose the whole eye. After head position
10 is okay, user can click the start button to start the training. When the user click the start button,
11 the following screen is shown.



- 12
13 In the upper left corner, it will have a counter indicate that the training will start. The user
14 should try to gaze the green point on the upper corner. Note that the green point will change

- 1 the position from upper left corner to upper right corner after certain gaze duration. After the
- 2 training complete, message will shown that training is completed.



- 3
- 4

6. Operation Phase

In this section, we will describe how to use the AppBlaster to control the electric appliances.

6.1. Move and selection

When the training is done, the following screen is shown



In the upper part of the screen, there is a camera view that contains the “ghost image” for the image capture during the training. The user should try to match the head with the “ghost image” for better gaze direction detection. In the upper left and right, there is a LED indicate that user’s gaze direction. When the user tries to gaze the left upper corner of the screen, the LED in the upper left corner should be on.



It is same for the upper right LED case.

In the lower part of the screen, it shows the current selected electric appliance. Note that for each electric appliance, there are two icons: the coloured icon to turning on the electric appliance and the black and white icon for turning off the electric appliance.


In order to trigger the current selected electric appliance, the user has to gaze the upper right corner continuously for a period of time. For change the selection, user has to gaze the left right corner for a period of time. Note that such period of time can be configured on settings.

User can also use finger touch to change the device by dragging from left to right or right to left on the lower part of the screen, and trigger the current selected item by long pressed the selected item.


6.2. Set the gaze duration

In the AppBlaster, the period of time required to gaze the corner continuously for triggering action can be configured. For longer such duration, there will be less false trigger, but it will also increase the difficulty and time for performing an action.

To change such duration, user can go to the settings screen.



Then user can select the gaze direction, the follow screen will be shown



Then user can click the time interval to select the time required for gaze on corner.




6.3. Lock screen

In the operation screen, besides the electric appliances, there is another icon for lock screen.



- 1
2 The lock screen can be used to lock the system until the user make correct gaze direction
3 sequence. When lock screen icon is selected, the following screen is shown.



- 4
5 In this lock screen, in order to unlock, user has to gaze right, left, right, right, left sequence.
6 When the gaze direction can be matched, the gaze direction icon will become gray.



- 7
8 When all the gaze direction is matched, the screen will be unlocked.